



☎ 55 157 506
90 157 560

🌐 www.gnie.com

✉ grtarek@yahoo.fr



Imm Emmeraude de Tunis Rue Mohamed Badra
Bureau A-2-8 Montplaisir, Tunis, 1073

2. Specific Operating Steps



Step 1: Disinfection

To ensure the hygiene and safety of the hair removal area, first clean and disinfect the hair removal area to avoid bacterial infection.



Step 2: Insertion

Use Insulated Needle
 • Insert the needle into the hair follicle at the correct angle.
 • Apply low-level current to destroy the hair root.



Step 3: Hot and Cold Hammer

• Cold Mode: Use before treatment to numb the skin and reduce pain.
 • Hot Mode: Use after treatment to soothe redness and irritation.



Step 4: Insulated Tweezers

• Gently remove the treated hair with insulated tweezers.
 • Ensure no residual hair remains in the follicle.



Step 5: Plasma Pen

• Use the plasma pen in low mode to stimulate collagen and promote healing.
 • Helps reduce redness and prevent ingrown hairs.

